

Snap-Off & SBCS

(Small Bone Compression Screw)



RESP  **NSE**
ORTH

A WishBone Medical Company

Response Ortho is a global
orthopaedic trauma solutions manufacturer
offering premium products
created under its founding principles of
innovation, excellence by design
and functional superiority





CONTENTS

- 3** Introductions
- 4** Snap-Off Screw
Product Features
- 5** Snap-Off Screw
Ordering Information
- 6** SBCS Product Features
- 7** SBCS Ordering Information

www.responseortho.com

Introduction


The Snap-Off and Small Bone Compression Screws can be used in combination with a drill, or the surgeon can choose to use the self drilling & self tapping feature of the screws.

The systems are indicated for fixation of small bone fractures, osteotomies and small joint fusion typically in the foot or hand.

These screw systems are designed to reduce the risk of the soft tissue irritation while providing accurate placement, controlled compression and stable fixation.

Indications

- Fixation of fractures, osteotomies small joint fusions typically within the small bones of the hand and the feet.



Implants

40-5082-11	2,0mm/11,0mm	4
40-5082-12	2,0mm/12,0mm	4
40-5082-13	2,0mm/13,0mm	4
40-5082-14	2,0mm/14,0mm	4



Instruments


00-1011-01	Quick Release Handle	1
------------	----------------------	---



00-2090-00	Snap - off Screw Driver	2
------------	-------------------------	---





00-8260-00	Snap - off Screw Tray	1
------------	-----------------------	---



Applications

- Hallux valgus
- Scaphoid Fractures
- Atkins Chevron Osteotomies (Smallbone Osteotomies)
- Fusion of MTP Joint
- Small Bone Fractures



Implants

40-5072-12	4mm/3mm x 12mm	3
40-5072-14	4mm/3mm x 14mm	3
40-5072-16	4mm/3mm x 16mm	4
40-5072-18	4mm/3mm x 18mm	4
40-5072-20	4mm/3mm x 20mm	4
40-5072-22	4mm/3mm x 22mm	3
40-5072-24	4mm/3mm x 24mm	3



Instruments

00-1011-01	Driver Handle Quick Release Handle	1
------------	------------------------------------	---



00-2084-00	Screw Holder	1
------------	--------------	---



00-0046-23	Drill Guide w / handle, 2.3mm	1
------------	-------------------------------	---



00-3223-01	2.3mm Can. Quick Release Drill	2
------------	--------------------------------	---



00-3339-01	3.3mm Can. Quick Release Counter Sink Short	2
------------	---	---



40-4101-15	Guide Wire, 1.0x150mm,Trocar	4
------------	------------------------------	---



00-0013-01	Wire Depth Gauge, Short	4
------------	-------------------------	---



00-2014-01	Quick Release Hex Driver Tip, 2.0mm, Short	2
------------	--	---




00-2014-01	Quick Release Hex Driver Tip, 2.0mm, Cannulated, Short	2
------------	--	---



00-0046-10	Drill Guide w / handle, for 1.0mm	1
------------	-----------------------------------	---



00-8270-00	ForeFoot Screw Tray	1
------------	---------------------	---



Snap-Off Screw SBCS

Small Bone Compression Screw